

Product datasheet for **TC400004**

Transient Overexpression Lysate of Myc-DDK-tGFP in HEK293T cells paraffin embedded cell block

Product data:

Product Type:	Cytosections
Description:	Transient Overexpression Lysate of Myc-DDK-tGFP in HEK293T cells paraffin embedded cell block
Control Type:	Negative
Species:	Human
Protein Source:	Synthetic
Tag:	Myc-DDK
Storage:	RT
Stability:	Stable for at least 1 year from receipt of products under proper storage and handling conditions.



[View online »](#)