

Product datasheet for **TF300647**

SUN1 Human shRNA Plasmid Kit (Locus ID 23353)

Product data:

Product Type:	shRNA Plasmids
Product Name:	SUN1 Human shRNA Plasmid Kit (Locus ID 23353)
Locus ID:	23353
Synonyms:	UNC84A
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	UNC84A - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 23353). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<u>NM_001130965</u> , <u>NM_001171944</u> , <u>NM_001171945</u> , <u>NM_001171946</u> , <u>NM_025154</u> , <u>NM_025154.1</u> , <u>NM_025154.2</u> , <u>NM_025154.3</u> , <u>NM_025154.4</u> , <u>NM_025154.5</u> , <u>NM_001130965.1</u> , <u>NM_001130965.2</u> , <u>NM_001171946.1</u> , <u>NM_001171945.1</u> , <u>NM_001171944.1</u> , <u>BC013613</u> , <u>BC113916</u> , <u>BC142707</u> , <u>BC151827</u> , <u>BC152418</u> , <u>BM689998</u> , <u>NM_001367633</u> , <u>NM_001367636</u> , <u>NM_001367638</u> , <u>NM_001367640</u> , <u>NM_001367643</u> , <u>NM_001367646</u> , <u>NM_001367655</u> , <u>NM_001367660</u> , <u>NM_001367666</u> , <u>NM_001367668</u> , <u>NM_001367670</u> , <u>NM_001367673</u> , <u>NM_001367675</u> , <u>NM_001367676</u> , <u>NM_001367678</u> , <u>NM_001367680</u> , <u>NM_001367681</u> , <u>NM_001367683</u> , <u>NM_001367684</u> , <u>NM_001367686</u> , <u>NM_001367690</u> , <u>NM_001367691</u> , <u>NM_001367692</u> , <u>NM_001367694</u> , <u>NM_001367696</u> , <u>NM_001367705</u> , <u>NM_001367708</u> , <u>NR_160282</u> , <u>NR_160283</u> , <u>NM_001367634</u> , <u>NM_001367635</u> , <u>NM_001367637</u> , <u>NM_001367639</u> , <u>NM_001367641</u> , <u>NM_001367642</u> , <u>NM_001367644</u> , <u>NM_001367645</u> , <u>NM_001367647</u> , <u>NM_001367648</u> , <u>NM_001367649</u> , <u>NM_001367651</u> , <u>NM_001367653</u> , <u>NM_001367658</u> , <u>NM_001367662</u> , <u>NM_001367664</u> , <u>NM_001367665</u> , <u>NM_001367667</u> , <u>NM_001367669</u> , <u>NM_001367671</u> , <u>NM_001367672</u> , <u>NM_001367674</u> , <u>NM_001367677</u> , <u>NM_001367679</u> , <u>NM_001367682</u> , <u>NM_001367685</u> , <u>NM_001367687</u> , <u>NM_001367688</u> , <u>NM_001367689</u> , <u>NM_001367693</u> , <u>NM_001367695</u> , <u>NM_001367697</u> , <u>NM_001367698</u> , <u>NM_001367699</u> , <u>NM_001367700</u> , <u>NM_001367701</u> , <u>NM_001367702</u> , <u>NM_001367703</u> , <u>NM_001367704</u> , <u>NM_001367706</u> , <u>NR_160281</u> , <u>NM_025154.6</u> , <u>NM_001171945.2</u> , <u>NM_001171944.2</u> , <u>NM_001130965.3</u> , <u>NM_001171946.2</u>



[View online »](#)

UniProt ID: [O94901](#)

Summary: This gene is a member of the unc-84 homolog family and encodes a nuclear envelope protein with an Unc84 (SUN) domain. The protein is involved in nuclear anchorage and migration. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jan 2019]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).