

## Product datasheet for **TF301508**

### SMARCA2 Human shRNA Plasmid Kit (Locus ID 6595)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SMARCA2 Human shRNA Plasmid Kit (Locus ID 6595)
Locus ID:	6595
Synonyms:	BAF190; BIS; BRM; hBRM; hSNF2a; NCBRS; SNF2; SNF2L2; SNF2LA; Sth1p; SWI2
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	SMARCA2 - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 6595). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">NM_001289396</a> , <a href="#">NM_001289397</a> , <a href="#">NM_001289398</a> , <a href="#">NM_001289399</a> , <a href="#">NM_001289400</a> , <a href="#">NM_003070</a> , <a href="#">NM_139045</a> , <a href="#">NM_139045.1</a> , <a href="#">NM_139045.2</a> , <a href="#">NM_139045.3</a> , <a href="#">NM_003070.1</a> , <a href="#">NM_003070.2</a> , <a href="#">NM_003070.3</a> , <a href="#">NM_003070.4</a> , <a href="#">NM_001289398.1</a> , <a href="#">NM_001289399.1</a> , <a href="#">NM_001289400.1</a> , <a href="#">NM_001289397.1</a> , <a href="#">NM_001289396.1</a> , <a href="#">BC040029</a> , <a href="#">BC068252</a> , <a href="#">BC156185</a> , <a href="#">BM921013</a> , <a href="#">NM_001289397.2</a> , <a href="#">NM_003070.5</a> , <a href="#">NM_139045.4</a> , <a href="#">NM_001289398.2</a>
UniProt ID:	<a href="#">P51531</a>
Summary:	The protein encoded by this gene is a member of the SWI/SNF family of proteins and is highly similar to the brahma protein of Drosophila. Members of this family have helicase and ATPase activities and are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SNF/SWI, which is required for transcriptional activation of genes normally repressed by chromatin. Alternatively spliced transcript variants encoding different isoforms have been found for this gene, which contains a trinucleotide repeat (CAG) length polymorphism. [provided by RefSeq, Jan 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).