

## Product datasheet for **TF305405**

### CHKL (CHKB) Human shRNA Plasmid Kit (Locus ID 1120)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CHKL (CHKB) Human shRNA Plasmid Kit (Locus ID 1120)
Locus ID:	1120
Synonyms:	CHETK; CHKL; CK; CKB; CKEKB; EK; EKB; MDCMC
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CHKB - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 1120). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">NM_005198</a> , <a href="#">NM_152253</a> , <a href="#">NM_005198.1</a> , <a href="#">NM_005198.2</a> , <a href="#">NM_005198.3</a> , <a href="#">NM_005198.4</a> , <a href="#">NM_152253.1</a> , <a href="#">BC082263</a> , <a href="#">BC024286</a> , <a href="#">BC037162</a> , <a href="#">BC050443</a> , <a href="#">BC056404</a> , <a href="#">BC101488</a> , <a href="#">BC113521</a> , <a href="#">NM_005198.5</a>
UniProt ID:	<a href="#">Q9Y259</a>
Summary:	Choline kinase (CK) and ethanolamine kinase (EK) catalyze the phosphorylation of choline/ethanolamine to phosphocholine/phosphoethanolamine. This is the first enzyme in the biosynthesis of phosphatidylcholine/phosphatidylethanolamine in all animal cells. The highly purified CKs from mammalian sources and their recombinant gene products have been shown to have EK activity also, indicating that both activities reside on the same protein. The choline kinase-like protein encoded by CHKL belongs to the choline/ethanolamine kinase family; however, its exact function is not known. Read-through transcripts are expressed from this locus that include exons from the downstream CPT1B locus. [provided by RefSeq, Jun 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).