

Product datasheet for TF308332

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

YAP1 Human shRNA Plasmid Kit (Locus ID 10413)

Product data:

Product Type: shRNA Plasmids

Product Name: YAP1 Human shRNA Plasmid Kit (Locus ID 10413)

Locus ID: 10413

Synonyms: COB1; YAP; YAP2; YAP65; YKI

Vector: pRFP-C-RS (TR30014)

E. coli Selection: Chloramphenicol (34 ug/ml)

Mammalian Cell

Selection:

Puromycin

Format: Retroviral plasmids

Components: YAP1 - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 10413).

5µg purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.

RefSeq: NM 001130145, NM 001195044, NM 001195045, NM 001282097, NM 001282098,

NM 001282099, NM 001282100, NM 001282101, NM 006106, NM 006106.1, NM 006106.2,

NM 006106.3, NM 006106.4, NM 001130145.1, NM 001130145.2, NM 001195045.1,

NM 001195044.1, NM 001282098.1, NM 001282097.1, NM 001282099.1, NM 001282100.1,

NM 001282101.1, BC038235, BC026212, NM 001195044.2, NM 001282097.2,

NM 001130145.3, NM 006106.5, NM 001195045.2

UniProt ID: P46937

Summary: This gene encodes a downstream nuclear effector of the Hippo signaling pathway which is

involved in development, growth, repair, and homeostasis. This gene is known to play a role in the development and progression of multiple cancers as a transcriptional regulator of this signaling pathway and may function as a potential target for cancer treatment. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by

RefSeq, Aug 2013]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u>.





Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).