

Product datasheet for **TF312078**

ITGB6 Human shRNA Plasmid Kit (Locus ID 3694)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ITGB6 Human shRNA Plasmid Kit (Locus ID 3694)
Locus ID:	3694
Synonyms:	AI1H
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ITGB6 - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 3694). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	NM_000888 , NM_001282353 , NM_001282354 , NM_001282355 , NM_001282388 , NM_001282389 , NM_001282390 , NM_000888.1 , NM_000888.2 , NM_000888.4 , NM_001282390.1 , NM_001282354.1 , NM_001282389.1 , NM_001282388.1 , NM_001282353.1 , BC121178 , BC143656 , NM_001282389.2 , NM_001282355.2 , NM_000888.5 , NM_001282388.2 , NM_001282353.2 , NM_001282390.2
UniProt ID:	P18564
Summary:	This gene encodes a protein that is a member of the integrin superfamily. Members of this family are adhesion receptors that function in signaling from the extracellular matrix to the cell. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. The encoded protein forms a dimer with an alpha v chain and this heterodimer can bind to ligands like fibronectin and transforming growth factor beta 1. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).