

## **Product datasheet for TF314721**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Apolipoprotein E (APOE) Human shRNA Plasmid Kit (Locus ID 348)**

**Product data:** 

**Product Type:** shRNA Plasmids

**Product Name:** Apolipoprotein E (APOE) Human shRNA Plasmid Kit (Locus ID 348)

Locus ID: 348

Synonyms: AD2; APO-E; ApoE4; LDLCQ5; LPG

**Vector:** pRFP-C-RS (TR30014)

E. coli Selection: Chloramphenicol (34 ug/ml)

**Mammalian Cell** 

Selection:

Puromycin

Format:

Retroviral plasmids

Components: APOE - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 348). 5µg

purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.

RefSeq: NM 000041, NM 001302688, NM 001302689, NM 001302690, NM 001302691, NM 000041.1,

NM 000041.2, NM 000041.3, NM 001302689.1, NM 001302690.1, NM 001302691.1, NM 001302688.1, BC003557, BC003557.1, BM042032, NM 001302688.2, NM 000041.4,

NM 001302689.2, NM 001302691.2

UniProt ID: P02649

**Summary:** The protein encoded by this gene is a major apoprotein of the chylomicron. It binds to a

specific liver and peripheral cell receptor, and is essential for the normal catabolism of triglyceride-rich lipoprotein constituents. This gene maps to chromosome 19 in a cluster with

the related apolipoprotein C1 and C2 genes. Mutations in this gene result in familial

dysbetalipoproteinemia, or type III hyperlipoproteinemia (HLP III), in which increased plasma cholesterol and triglycerides are the consequence of impaired clearance of chylomicron and

VLDL remnants. [provided by RefSeq, Jun 2016]

**shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our custom shRNA service.







## Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).