

## Product datasheet for **TF317099**

### USP21 Human shRNA Plasmid Kit (Locus ID 27005)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	USP21 Human shRNA Plasmid Kit (Locus ID 27005)
Locus ID:	27005
Synonyms:	USP16; USP23
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	USP21 - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 27005). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">NM_001014443</a> , <a href="#">NM_001319848</a> , <a href="#">NM_012475</a> , <a href="#">NM_016572</a> , <a href="#">NM_001319847</a> , <a href="#">NM_012475.1</a> , <a href="#">NM_012475.2</a> , <a href="#">NM_012475.3</a> , <a href="#">NM_001014443.1</a> , <a href="#">NM_001014443.2</a> , <a href="#">BC090946</a> , <a href="#">BC090946.1</a> , <a href="#">BC003130</a> , <a href="#">BC136291</a> , <a href="#">BM742944</a> , <a href="#">NM_001014443.3</a> , <a href="#">NM_012475.5</a>
UniProt ID:	<a href="#">Q9UK80</a>
Summary:	This gene encodes a member of the C19 peptidase family, also known as family 2 of ubiquitin carboxy-terminal hydrolases. The encoded protein cleaves ubiquitin from ubiquitinated proteins for recycling in intracellular protein degradation. The encoded protein is also able to release NEDD8, a ubiquitin-like protein, from NEDD8-conjugated proteins. This gene has been referred to as USP16 and USP23 but is now known as USP21. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).