

Product datasheet for TF320279

BubR1 (BUB1B) Human shRNA Plasmid Kit (Locus ID 701)

Product data:

Product Type: shRNA Plasmids **Product Name:** BubR1 (BUB1B) Human shRNA Plasmid Kit (Locus ID 701) Locus ID: 701 Bub1A; BUB1beta; BUBR1; hBUBR1; MAD3L; MVA1; SSK1 Synonyms: pRFP-C-RS (TR30014) Vector: E. coli Selection: Chloramphenicol (34 ug/ml) Mammalian Cell Puromycin Selection: Format: **Retroviral plasmids Components:** BUB1B - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 701). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free. NM 001211, NM 001211.1, NM 001211.2, NM 001211.3, NM 001211.4, NM 001211.5, RefSeq: BC018739, BC018739.1 **UniProt ID:** 060566 Summary: This gene encodes a kinase involved in spindle checkpoint function. The protein has been localized to the kinetochore and plays a role in the inhibition of the anaphase-promoting complex/cyclosome (APC/C), delaying the onset of anaphase and ensuring proper chromosome segregation. Impaired spindle checkpoint function has been found in many forms of cancer. [provided by RefSeq, Jul 2008] These shRNA constructs were designed against multiple splice variants at this gene locus. To shRNA Design: be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our custom shRNA service.



View online 🤉

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GRIGENE BubR1 (BUB1B) Human shRNA Plasmid Kit (Locus ID 701) – TF320279

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US