

## Product datasheet for **TF320285**

### **CAMK2G Human shRNA Plasmid Kit (Locus ID 818)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	CAMK2G Human shRNA Plasmid Kit (Locus ID 818)
Locus ID:	818
Synonyms:	CAMK; CAMK-II; CAMKG; MRD59
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CAMK2G - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 818). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">NM_001204492</a> , <a href="#">NM_001222</a> , <a href="#">NM_001320898</a> , <a href="#">NM_172169</a> , <a href="#">NM_172170</a> , <a href="#">NM_172171</a> , <a href="#">NM_172172</a> , <a href="#">NM_172173</a> , <a href="#">NM_172169.1</a> , <a href="#">NM_172169.2</a> , <a href="#">NM_172173.1</a> , <a href="#">NM_172173.2</a> , <a href="#">NM_001222.1</a> , <a href="#">NM_001222.2</a> , <a href="#">NM_001222.3</a> , <a href="#">NM_172171.1</a> , <a href="#">NM_172171.2</a> , <a href="#">NM_172170.1</a> , <a href="#">NM_172170.2</a> , <a href="#">NM_172170.3</a> , <a href="#">NM_001204492.1</a> , <a href="#">NM_172172.1</a> , <a href="#">BC034044</a> , <a href="#">BC012795</a> , <a href="#">BC021269</a> , <a href="#">BC037928</a> , <a href="#">NM_001367547</a> , <a href="#">NM_001367516</a> , <a href="#">NM_001367518</a> , <a href="#">NM_001367525</a> , <a href="#">NM_001367526</a> , <a href="#">NM_001367529</a> , <a href="#">NM_001367534</a> , <a href="#">NM_001367537</a> , <a href="#">NM_001367539</a> , <a href="#">NM_001367540</a> , <a href="#">NM_001367542</a> , <a href="#">NM_001367544</a> , <a href="#">NM_001367545</a> , <a href="#">NR_160042</a> , <a href="#">NR_160045</a> , <a href="#">NR_160263</a> , <a href="#">NM_001367514</a> , <a href="#">NM_001367517</a> , <a href="#">NM_001367519</a> , <a href="#">NM_001367520</a> , <a href="#">NM_001367521</a> , <a href="#">NM_001367522</a> , <a href="#">NM_001367523</a> , <a href="#">NM_001367524</a> , <a href="#">NM_001367527</a> , <a href="#">NM_001367528</a> , <a href="#">NM_001367530</a> , <a href="#">NM_001367531</a> , <a href="#">NM_001367532</a> , <a href="#">NM_001367533</a> , <a href="#">NM_001367535</a> , <a href="#">NM_001367536</a> , <a href="#">NM_001367538</a> , <a href="#">NM_001367541</a> , <a href="#">NM_001367543</a> , <a href="#">NM_001367546</a> , <a href="#">NM_001367548</a> , <a href="#">NR_160040</a> , <a href="#">NR_160041</a> , <a href="#">NR_160044</a> , <a href="#">NR_160046</a> , <a href="#">NR_160047</a> , <a href="#">NM_001222.4</a> , <a href="#">NM_172173.3</a> , <a href="#">NM_172169.3</a> , <a href="#">NM_172170.5</a>
UniProt ID:	<a href="#">Q13555</a>



[View online »](#)

<b>Summary:</b>	<p>The product of this gene is one of the four subunits of an enzyme which belongs to the serine/threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. Calcium signaling is crucial for several aspects of plasticity at glutamatergic synapses. In mammalian cells the enzyme is composed of four different chains: alpha, beta, gamma, and delta. The product of this gene is a gamma chain. Many alternatively spliced transcripts encoding different isoforms have been described but the full-length nature of all the variants has not been determined.[provided by RefSeq, Mar 2011]</p>
<b>shRNA Design:</b>	<p>These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a>.</p>
<b>Performance Guaranteed:</b>	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>