

## Product datasheet for **TF320334**

### Eph receptor A7 (EPHA7) Human shRNA Plasmid Kit (Locus ID 2045)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Eph receptor A7 (EPHA7) Human shRNA Plasmid Kit (Locus ID 2045)
Locus ID:	2045
Synonyms:	EHK-3; EHK3; EK11; HEK11
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	EPHA7 - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 2045). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">NM_001288629</a> , <a href="#">NM_001288630</a> , <a href="#">NM_004440</a> , <a href="#">NM_004440.1</a> , <a href="#">NM_004440.2</a> , <a href="#">NM_004440.3</a> , <a href="#">NM_001288629.1</a> , <a href="#">BC126125</a> , <a href="#">BC027940</a> , <a href="#">BC126151</a> , <a href="#">BC143857</a> , <a href="#">BC143858</a> , <a href="#">BM677016</a> , <a href="#">NM_001288629.2</a> , <a href="#">NM_001288630.2</a> , <a href="#">NM_004440.4</a>
UniProt ID:	<a href="#">Q15375</a>
Summary:	This gene belongs to the ephrin receptor subfamily of the protein-tyrosine kinase family. EPH and EPH-related receptors have been implicated in mediating developmental events, particularly in the nervous system. Receptors in the EPH subfamily typically have a single kinase domain and an extracellular region containing a Cys-rich domain and 2 fibronectin type III repeats. The ephrin receptors are divided into 2 groups based on the similarity of their extracellular domain sequences and their affinities for binding ephrin-A and ephrin-B ligands. Increased expression of this gene is associated with multiple forms of carcinoma. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).