

Product datasheet for **TF320382**

HNF 4 alpha (HNF4A) Human shRNA Plasmid Kit (Locus ID 3172)

Product data:

Product Type:	shRNA Plasmids
Product Name:	HNF 4 alpha (HNF4A) Human shRNA Plasmid Kit (Locus ID 3172)
Locus ID:	3172
Synonyms:	FRTS4; HNF4; HNF4a7; HNF4a8; HNF4a9; HNF4alpha; MODY; MODY1; NR2A1; NR2A21; TCF; TCF-14; TCF14
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	HNF4A - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 3172). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	NM_000457 , NM_001030003 , NM_001030004 , NM_001258355 , NM_001287182 , NM_001287183 , NM_001287184 , NM_175914 , NM_178849 , NM_178850 , NM_001030003.1 , NM_001030003.2 , NM_178849.1 , NM_178849.2 , NM_175914.3 , NM_175914.4 , NM_001030004.1 , NM_001030004.2 , NM_000457.1 , NM_000457.2 , NM_000457.3 , NM_000457.4 , NM_178850.1 , NM_178850.2 , NM_001258355.1 , NM_001287184.1 , NM_001287182.1 , NM_001287183.1 , BC137539 , BC137540 , NM_000457.5 , NM_178850.3 , NM_001030003.3 , NM_001287182.2 , NM_178849.3 , NM_001287184.2 , NM_001258355.2 , NM_001030004.3
UniProt ID:	P41235
Summary:	The protein encoded by this gene is a nuclear transcription factor which binds DNA as a homodimer. The encoded protein controls the expression of several genes, including hepatocyte nuclear factor 1 alpha, a transcription factor which regulates the expression of several hepatic genes. This gene may play a role in development of the liver, kidney, and intestines. Mutations in this gene have been associated with monogenic autosomal dominant non-insulin-dependent diabetes mellitus type I. Alternative splicing of this gene results in multiple transcript variants encoding several different isoforms. [provided by RefSeq, Apr 2012]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).