

## Product datasheet for **TF320537**

### **MST1 (STK4) Human shRNA Plasmid Kit (Locus ID 6789)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	MST1 (STK4) Human shRNA Plasmid Kit (Locus ID 6789)
Locus ID:	6789
Synonyms:	KRS2; MST1; YSK3
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	STK4 - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 6789). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">BC005231</a> , <a href="#">NM_006282</a> , <a href="#">NM_001352385</a> , <a href="#">NR_147974</a> , <a href="#">NR_147975</a> , <a href="#">NM_006282.1</a> , <a href="#">NM_006282.2</a> , <a href="#">NM_006282.3</a> , <a href="#">BC005231.1</a> , <a href="#">BC093768</a> , <a href="#">BC093768.1</a> , <a href="#">BC015332</a> , <a href="#">BC029511</a> , <a href="#">BC038850</a> , <a href="#">BC039023</a> , <a href="#">BC042861</a> , <a href="#">BC058916</a> , <a href="#">BC065215</a> , <a href="#">BM450795</a> , <a href="#">BM455776</a> , <a href="#">BM459914</a> , <a href="#">NM_006282.5</a>
UniProt ID:	<a href="#">Q13043</a>
Summary:	The protein encoded by this gene is a cytoplasmic kinase that is structurally similar to the yeast Ste20p kinase, which acts upstream of the stress-induced mitogen-activated protein kinase cascade. The encoded protein can phosphorylate myelin basic protein and undergoes autophosphorylation. A caspase-cleaved fragment of the encoded protein has been shown to be capable of phosphorylating histone H2B. The particular phosphorylation catalyzed by this protein has been correlated with apoptosis, and it's possible that this protein induces the chromatin condensation observed in this process. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).