

## Product datasheet for **TF320553**

### TLR2 Human shRNA Plasmid Kit (Locus ID 7097)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	TLR2 Human shRNA Plasmid Kit (Locus ID 7097)
Locus ID:	7097
Synonyms:	CD282; TIL4
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	TLR2 - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 7097). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">NM_001318787</a> , <a href="#">NM_001318789</a> , <a href="#">NM_001318790</a> , <a href="#">NM_001318791</a> , <a href="#">NM_001318793</a> , <a href="#">NM_001318795</a> , <a href="#">NM_001318796</a> , <a href="#">NM_003264</a> , <a href="#">NM_003264.1</a> , <a href="#">NM_003264.2</a> , <a href="#">NM_003264.3</a> , <a href="#">NM_003264.4</a> , <a href="#">BC033756</a> , <a href="#">BC033756.1</a> , <a href="#">NM_003264.5</a>
UniProt ID:	<a href="#">O60603</a>
Summary:	The protein encoded by this gene is a member of the Toll-like receptor (TLR) family which plays a fundamental role in pathogen recognition and activation of innate immunity. TLRs are highly conserved from Drosophila to humans and share structural and functional similarities. This protein is a cell-surface protein that can form heterodimers with other TLR family members to recognize conserved molecules derived from microorganisms known as pathogen-associated molecular patterns (PAMPs). Activation of TLRs by PAMPs leads to an up-regulation of signaling pathways to modulate the host's inflammatory response. This gene is also thought to promote apoptosis in response to bacterial lipoproteins. This gene has been implicated in the pathogenesis of several autoimmune diseases. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).