

Product datasheet for **TF320665**

TNIK Human shRNA Plasmid Kit (Locus ID 23043)

Product data:

Product Type:	shRNA Plasmids
Product Name:	TNIK Human shRNA Plasmid Kit (Locus ID 23043)
Locus ID:	23043
Synonyms:	MRT54
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	KIAA0551 - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 23043). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	NM_001161560 , NM_001161561 , NM_001161562 , NM_001161563 , NM_001161564 , NM_001161565 , NM_001161566 , NM_015028 , NR_027767 , NM_015028.1 , NM_015028.2 , NM_015028.3 , NM_001161566.1 , NM_001161566.2 , NM_001161565.1 , NM_001161565.2 , NM_001161564.1 , NM_001161564.2 , NM_001161563.1 , NM_001161563.2 , NM_001161562.1 , NM_001161562.2 , NM_001161561.1 , NM_001161561.2 , NM_001161560.1 , NM_001161560.2 , BC018093 , BC019081 , BC055427 , BC150256
UniProt ID:	Q9UKE5
Summary:	Wnt signaling plays important roles in carcinogenesis and embryonic development. The protein encoded by this gene is a serine/threonine kinase that functions as an activator of the Wnt signaling pathway. Mutations in this gene are associated with an autosomal recessive form of cognitive disability. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).