

## Product datasheet for **TF320698**

### CrkRS (CDK12) Human shRNA Plasmid Kit (Locus ID 51755)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CrkRS (CDK12) Human shRNA Plasmid Kit (Locus ID 51755)
Locus ID:	51755
Synonyms:	CRK7; CRKR; CRKRS
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CDK12 - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 51755). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">NM_015083</a> , <a href="#">NM_016507</a> , <a href="#">NM_016507.1</a> , <a href="#">NM_016507.2</a> , <a href="#">NM_016507.3</a> , <a href="#">NM_015083.1</a> , <a href="#">NM_015083.2</a> , <a href="#">BC140854</a> , <a href="#">BC150265</a> , <a href="#">BM557102</a> , <a href="#">NM_015083.3</a> , <a href="#">NM_016507.4</a>
UniProt ID:	<a href="#">Q9NYV4</a>
Summary:	Cyclin-dependent kinase that phosphorylates the C-terminal domain (CTD) of the large subunit of RNA polymerase II (POLR2A), thereby acting as a key regulator of transcription elongation. Regulates the expression of genes involved in DNA repair and is required for the maintenance of genomic stability. Preferentially phosphorylates 'Ser-5' in CTD repeats that are already phosphorylated at 'Ser-7', but can also phosphorylate 'Ser-2'. Required for RNA splicing, possibly by phosphorylating SRSF1/SF2. Involved in regulation of MAP kinase activity, possibly leading to affect the response to estrogen inhibitors.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).