

## Product datasheet for **TF320741**

### SGK196 (POMK) Human shRNA Plasmid Kit (Locus ID 84197)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | shRNA Plasmids   |
| Product Name:             | SGK196 (POMK) Human shRNA Plasmid Kit (Locus ID 84197)   |
| Locus ID:                 | 84197  |
| Synonyms:                 | MDDGA12; MDDGC12; SGK196   |
| Vector:                   | pRFP-C-RS (TR30014)  |
| E. coli Selection:        | Chloramphenicol (34 ug/ml)   |
| Mammalian Cell Selection: | Puromycin  |
| Format:                   | Retroviral plasmids  |
| Components:               | SGK196 - Human, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 84197). 5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.  |
| RefSeq:                   | <a href="#">NM_001277971</a> , <a href="#">NM_032237</a> , <a href="#">NM_032237.1</a> , <a href="#">NM_032237.2</a> , <a href="#">NM_032237.3</a> , <a href="#">NM_032237.4</a> , <a href="#">NM_001277971.1</a> , <a href="#">BC113703</a> , <a href="#">BC101548</a> , <a href="#">NM_032237.5</a>  |
| UniProt ID:               | <a href="#">Q9H5K3</a>   |
| Summary:                  | This gene encodes a protein that may be involved in the presentation of the laminin-binding O-linked carbohydrate chain of alpha-dystroglycan (α-DG), which forms transmembrane linkages between the extracellular matrix and the exoskeleton. Some pathogens use this O-linked carbohydrate unit for host entry. Loss of function compound heterozygous mutations in this gene were found in a human patient affected by the Walker-Warburg syndrome (WWS) phenotype. Mice lacking this gene contain misplaced neurons (heterotopia) in some regions of the brain, possibly from defects in neuronal migration. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, May 2013] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .   |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).