

Product datasheet for **TF500218**

Bmpr2 Mouse shRNA Plasmid (Locus ID 12168)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Bmpr2 Mouse shRNA Plasmid (Locus ID 12168)
Locus ID:	12168
Synonyms:	2610024H22Rik; AL117858; AW546137; BB189135; BM; BMP; BMP-; BMP-2; BMPR-2; BMPR-II; BMPRII; BRK-3; Gm20272
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Bmpr2 - Mouse, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 12168). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	NM_007561 , NM_007561.1 , NM_007561.2 , NM_007561.3 , NM_007561.4 , BC148339 , BC156868 , BM940059
UniProt ID:	Q35607
Summary:	This gene encodes a serine/threonine kinase that functions as a receptor for bone morphogenetic proteins (BMPs). The encoded protein is a type II receptor that binds extracellular BMPs and forms a complex of two type II and two type I receptors at the cell membrane. This complex signals downstream to activate SMAD transcriptional regulators. This signaling is important during embryonic development. Mutations in this gene can cause pulmonary hypertension. [provided by RefSeq, Mar 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).