

Product datasheet for **TF500244**

Cacna1e Mouse shRNA Plasmid (Locus ID 12290)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Cacna1e Mouse shRNA Plasmid (Locus ID 12290)
Locus ID:	12290
Synonyms:	A430040I15; alp; alpha1E; BII; Cach6; Cacn1a6; Cav2; Cav2.3; Cchr; Cchra1
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Cacna1e - Mouse, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 12290). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	NM_009782 , NM_009782.1 , NM_009782.2 , NM_009782.3 , BC150750
Summary:	This gene encodes an integral membrane protein that belongs to the calcium channel alpha-1 subunits family. Voltage-sensitive calcium channels mediate the entry of calcium ions into excitable cells and are also involved in a variety of calcium-dependent processes. Voltage-dependent calcium channels are multi-subunit complexes, comprised of alpha-1, alpha-2, beta and delta subunits in a 1:1:1:1 ratio. The isoform alpha-1E gives rise to R-type calcium currents and belongs to the high-voltage activated group. Calcium channels containing the alpha-1E subunit may be involved in the modulation of neuronal firing patterns, an important component of information processing. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).