

## Product datasheet for **TF500797**

### Gfra1 Mouse shRNA Plasmid (Locus ID 14585)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Gfra1 Mouse shRNA Plasmid (Locus ID 14585)
Locus ID:	14585
Synonyms:	AU042498; GDNFR; GFRalpha-1
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Gfra1 - Mouse, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 14585). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">BC054378</a> , <a href="#">NM_001285457</a> , <a href="#">NM_010279</a> , <a href="#">NR_104342</a> , <a href="#">NM_010279.1</a> , <a href="#">NM_010279.2</a> , <a href="#">NM_010279.3</a> , <a href="#">NM_001285457.1</a> , <a href="#">NM_001285457.2</a> , <a href="#">BC011532</a> , <a href="#">BC040251</a> , <a href="#">BC044783</a>
UniProt ID:	<a href="#">P97785</a>
Summary:	This gene encodes a transmembrane protein that functions as the receptor for glial cell line derived neurotrophic factor (GDNF). The encoded protein undergoes proteolytic processing to generate a glycosylphosphatidylinositol-anchored cell surface coreceptor that forms a complex with the Ret tyrosine kinase in GDNF signaling pathway. Mice lacking the encoded protein exhibit deficits in the kidneys, the enteric nervous system, and spinal motor and sensory neurons similar mice deficient in GDNF or Ret. [provided by RefSeq, Jul 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).