

## Product datasheet for **TF509551**

### Ifnar1 Mouse shRNA Plasmid (Locus ID 15975)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Ifnar1 Mouse shRNA Plasmid (Locus ID 15975)
Locus ID:	15975
Synonyms:	CD118; Ifar; Ifnar; Ifrc; Infar
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Ifnar1 - Mouse, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 15975). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">BC043935</a> , <a href="#">BC052217</a> , <a href="#">BC052429</a> , <a href="#">NM_010508</a> , <a href="#">NM_010508.1</a> , <a href="#">NM_010508.2</a>
UniProt ID:	<a href="#">P33896</a>
Summary:	Component of the receptor for type I interferons, including interferons alpha, IFNB1 and IFNW1 (PubMed:1533935, PubMed:14532120, PubMed:23872679). Functions in general as heterodimer with IFNAR2 (By similarity). Type I interferon binding activates the JAK-STAT signaling cascade, and triggers tyrosine phosphorylation of a number of proteins including JAKs, TYK2, STAT proteins and the IFNR alpha- and beta-subunits themselves (PubMed:14532120). Can form an active IFNB1 receptor by itself and activate a signaling cascade that does not involve activation of the JAK-STAT pathway (PubMed:23872679). Contributes to modulate the innate immune response to bacterial lipopolysaccharide (PubMed:23872679).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).