

## Product datasheet for **TF511770**

### Spire2 Mouse shRNA Plasmid (Locus ID 234857)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Spire2 Mouse shRNA Plasmid (Locus ID 234857)
Locus ID:	234857
Synonyms:	BC026502; Spir-2; Spir2
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	SPIRE2 - Mouse, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 234857). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">BC026502</a> , <a href="#">BC049152</a> , <a href="#">NM_172287</a> , <a href="#">NM_172287.1</a> , <a href="#">NM_172287.2</a> , <a href="#">BC026502.1</a> , <a href="#">BC054072</a>
UniProt ID:	<a href="#">Q8K1S6</a>
Summary:	Acts as an actin nucleation factor, remains associated with the slow-growing pointed end of the new filament (PubMed:21620703, PubMed:21983562). Involved in intracellular vesicle transport along actin fibers, providing a novel link between actin cytoskeleton dynamics and intracellular transport (PubMed:21983562). Required for asymmetric spindle positioning and asymmetric cell division during oocyte meiosis (PubMed:21620703). Required for normal formation of the cleavage furrow and for polar body extrusion during female germ cell meiosis (PubMed:21620703). Also acts in the nucleus: together with SPIRE1 and SPIRE2, promotes assembly of nuclear actin filaments in response to DNA damage in order to facilitate movement of chromatin and repair factors after DNA damage (By similarity). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).