

## Product datasheet for **TF512980**

### Nat6 Mouse shRNA Plasmid (Locus ID 56441)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Nat6 Mouse shRNA Plasmid (Locus ID 56441)
Locus ID:	56441
Synonyms:	AI225910; Fus2
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Nat6 - Mouse, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 56441). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	<a href="#">BC026545</a> , <a href="#">NM_019750</a> , <a href="#">NM_019750.1</a> , <a href="#">NM_019750.2</a> , <a href="#">NM_019750.3</a> , <a href="#">NM_019750.4</a>
UniProt ID:	<a href="#">Q9R123</a>
Summary:	N-alpha-acetyltransferase that specifically mediates the acetylation of the acidic amino terminus of processed forms of beta- and gamma-actin (ACTB and ACTG, respectively). N-terminal acetylation of processed beta- and gamma-actin regulates actin filament depolymerization and elongation. In vivo, preferentially displays N-terminal acetyltransferase activity towards acid N-terminal sequences starting with Asp-Asp-Asp and Glu-Glu-Glu. In vitro, shows high activity towards Met-Asp-Glu-Leu and Met-Asp-Asp-Asp. May act as a tumor suppressor.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).