

Product datasheet for **TF515402**

Ddx17 Mouse shRNA Plasmid (Locus ID 67040)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Ddx17 Mouse shRNA Plasmid (Locus ID 67040)
Locus ID:	67040
Synonyms:	2610007K22Rik; A430025E01Rik; AI047725; C80929; Gm926; p7; p72
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Ddx17 - Mouse, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 67040). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	BC062910 , BC096036 , NM_001040187 , NM_152806 , NM_199079 , NM_199080 , NM_199079.1 , NM_001040187.1 , NM_199080.1 , NM_199080.2 , NM_152806.1 , NM_152806.2 , NM_152806.3 , BC012249 , BC027758 , BC031802 , BC038378
UniProt ID:	Q501J6
Summary:	This gene encodes the mouse homolog of human DEAD box polypeptide 17. DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD). RNA helicases of the DEAD-box family are involved in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Alternative splicing of this gene results in several transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).