

Product datasheet for **TF515799**

Gla3 Mouse shRNA Plasmid (Locus ID 110304)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Gla3 Mouse shRNA Plasmid (Locus ID 110304)
Locus ID:	110304
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Gla3 - Mouse, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 110304). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	NM_080438 , NM_080438.1 , NM_080438.2 , BC140461 , BC148711 , NM_001368774
UniProt ID:	Q91XP5
Summary:	Glycine receptors are ligand-gated chloride channels. Channel opening is triggered by extracellular glycine (PubMed:15131310, PubMed:20978350). Channel characteristics depend on the subunit composition; heteropentameric channels display faster channel closure (By similarity). Plays an important role in the down-regulation of neuronal excitability. Contributes to the generation of inhibitory postsynaptic currents (PubMed:15131310). Contributes to increased pain perception in response to increased prostaglandin E2 levels (PubMed:15131310). Plays a role in the regulation of breathing rhythm, especially of the duration of the postinspiratory phase (PubMed:20978350). Plays a role in cellular responses to ethanol (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).