

Product datasheet for **TF517150**

Tfcp2l1 Mouse shRNA Plasmid (Locus ID 81879)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Tfcp2l1 Mouse shRNA Plasmid (Locus ID 81879)
Locus ID:	81879
Synonyms:	1810030F05Rik; 4932442M07Rik; AA575098; Cp2l1; Crtr-1; D930018N21Rik; LBP-9; Tfcp2l1
Vector:	pRFP-C-RS (TR30014)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Tfcp2l1 - Mouse, 4 unique 29mer shRNA constructs in retroviral RFP vector(Gene ID = 81879). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRFP-C-RS Vector, TR30015, included for free.
RefSeq:	BC116663 , NM_023755 , NM_023755.1 , NM_023755.2
UniProt ID:	Q3UNW5
Summary:	Transcription factor that facilitates establishment and maintenance of pluripotency in embryonic stem cells (ESCs) (PubMed:23942233, PubMed:26321140). With Klf2, acts as the major effector of self-renewal that mediates induction of pluripotency downstream of LIF/Stat3 and Wnt/beta-catenin signaling (PubMed:23942238, PubMed:23942233, PubMed:26321140). Required for normal duct development in the salivary gland and kidney (PubMed:17079272). Coordinates the development of the kidney collecting ducts intercalated (IC) and principal (PC) cells, which regulate acid-base and salt-water homeostasis, respectively (PubMed:28577314). Regulates the expression of IC genes including subunits B1 and D2 of the V-ATPase complex, Oxgr1, Ca12, Slc4a1, Aqp6 and IC-specific transcription factor Foxi1 (PubMed:28577314). Regulates also the expression of Jag1 and subsequent notch signaling in the collecting duct (PubMed:28577314). Jag1 initiates notch signaling in PCs but inhibits notch signaling in ICs (PubMed:28577314). Acts as a transcriptional suppressor that may suppress UBP1-mediated transcriptional activation (PubMed:11073954). Modulates the placental expression of CYP11A1 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).