

Product datasheet for TF81001

TurboFectin 8.0 Transfection Reagent (1 mL in 1 vial)

Product data:

Product Type: Transfection Reagents

Description: TurboFectin Transfection Reagent (1 mL in 1 vial)


Features:

- Wide cell spectrum: validated on >100 cell lines and primary cells, [view list](#)
- Very low cell toxicity
- Simplified use: compatible with serum and antibiotics in the media

Storage: Upon receipt, please store the TurboFectin Transfection Reagent at +4°C.

Stability: The transfection reagent is stable for one (1) year from date of shipping when stored at +4°C.

Product images:



[View online »](#)