


Product datasheet for **TF81005**

TurboFectin 8.0 Transfection Reagent (5 mL in 5 vials)

Product data:

- Product Type:** Transfection Reagents
- Description:** TurboFectin Transfection Reagent (5 mL in 5 vials)
- Features:**
- Wide cell spectrum: validated on >100 cell lines and primary cells, [view list](#)
 - Very low cell toxicity
 - Simplified use: compatible with serum and antibiotics in the media
- Storage:** Upon receipt, please store the TurboFectin Transfection Reagent at +4°C.
- Stability:** The transfection reagent is stable for one (1) year from date of shipping when stored at +4°C.

Product images:



[View online »](#)