

## Product datasheet for **TG301162**

### TDRD5 Human shRNA Plasmid Kit (Locus ID 163589)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	TDRD5 Human shRNA Plasmid Kit (Locus ID 163589)
Locus ID:	163589
Synonyms:	TUDOR3
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	TDRD5 - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 163589). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	<a href="#">NM_001199085</a> , <a href="#">NM_001199089</a> , <a href="#">NM_001199091</a> , <a href="#">NM_001199092</a> , <a href="#">NM_173533</a> , <a href="#">NM_001349923</a> , <a href="#">NM_173533.1</a> , <a href="#">NM_173533.2</a> , <a href="#">NM_173533.3</a> , <a href="#">NM_001199085.1</a> , <a href="#">NM_001199085.2</a> , <a href="#">NM_001199092.1</a> , <a href="#">NM_001199091.1</a> , <a href="#">NM_001199089.1</a> , <a href="#">NM_001199089.2</a> , <a href="#">BC144058</a> , <a href="#">BC043368</a> , <a href="#">BC065707</a> , <a href="#">BC089387</a> , <a href="#">BC130532</a> , <a href="#">NM_001199089.3</a> , <a href="#">NM_001199085.3</a>
UniProt ID:	<a href="#">Q8NAT2</a>
Summary:	Required during spermiogenesis to participate in the repression transposable elements and prevent their mobilization, which is essential for the germline integrity. Probably acts via the piRNA metabolic process, which mediates the repression of transposable elements during meiosis by forming complexes composed of piRNAs and Piwi proteins and govern the methylation and subsequent repression of transposons. Required for chromatoid body (CB) assembly (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).