

## Product datasheet for **TG304713**

### EZH2 Human shRNA Plasmid Kit (Locus ID 2146)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	EZH2 Human shRNA Plasmid Kit (Locus ID 2146)
Locus ID:	2146
Synonyms:	ENX-1; ENX1; EZH2b; KMT6; KMT6A; WVS; WVS2
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	EZH2 - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 2146). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	<a href="#">NM_001203247</a> , <a href="#">NM_001203248</a> , <a href="#">NM_001203249</a> , <a href="#">NM_004456</a> , <a href="#">NM_152998</a> , <a href="#">NM_004456.1</a> , <a href="#">NM_004456.2</a> , <a href="#">NM_004456.3</a> , <a href="#">NM_004456.4</a> , <a href="#">NM_152998.1</a> , <a href="#">NM_152998.2</a> , <a href="#">NM_001203249.1</a> , <a href="#">NM_001203248.1</a> , <a href="#">NM_001203247.1</a> , <a href="#">BC010858</a> , <a href="#">BC010858.2</a> , <a href="#">BM974565</a> , <a href="#">NM_001203249.2</a> , <a href="#">NM_001203247.2</a> , <a href="#">NM_004456.5</a> , <a href="#">NM_001203248.2</a> , <a href="#">NM_152998.3</a>
UniProt ID:	<a href="#">Q15910</a>
Summary:	This gene encodes a member of the Polycomb-group (PcG) family. PcG family members form multimeric protein complexes, which are involved in maintaining the transcriptional repressive state of genes over successive cell generations. This protein associates with the embryonic ectoderm development protein, the VAV1 oncoprotein, and the X-linked nuclear protein. This protein may play a role in the hematopoietic and central nervous systems. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Feb 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).