

## Product datasheet for **TG304935**

### **DOCK5 Human shRNA Plasmid Kit (Locus ID 80005)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	DOCK5 Human shRNA Plasmid Kit (Locus ID 80005)
<b>Locus ID:</b>	80005
<b>Vector:</b>	pGFP-V-RS (TR30007)
<b>E. coli Selection:</b>	Kanamycin
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Retroviral plasmids
<b>Components:</b>	DOCK5 - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 80005). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
<b>RefSeq:</b>	<a href="#">BC011877</a> , <a href="#">NM_001322810</a> , <a href="#">NM_024940</a> , <a href="#">NM_024940.2</a> , <a href="#">NM_024940.3</a> , <a href="#">NM_024940.4</a> , <a href="#">NM_024940.5</a> , <a href="#">NM_024940.6</a> , <a href="#">NM_024940.7</a> , <a href="#">BC015352</a> , <a href="#">BC041761</a> , <a href="#">BC057384</a> , <a href="#">BC137175</a> , <a href="#">BC137176</a> , <a href="#">BM450222</a>
<b>UniProt ID:</b>	<a href="#">Q9H7D0</a>
<b>Summary:</b>	This gene encodes a member of the dedicator of cytokinesis protein family. Members of this family act as guanine nucleotide exchange factors for small Rho family G proteins. The protein encoded by this gene is thought to associate with adaptors CRK and CRKL, and function in regulation of intestinal epithelial cell spreading and migration on collagen IV. Similar proteins in mouse and zebrafish also function in myoblast fusion. [provided by RefSeq, Oct 2016]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).