

## Product datasheet for **TG308825**

### **TIAM1 Human shRNA Plasmid Kit (Locus ID 7074)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	TIAM1 Human shRNA Plasmid Kit (Locus ID 7074)
Locus ID:	7074
Synonyms:	TIAM-1
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	TIAM1 - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 7074). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	<a href="#">NM_003253</a> , <a href="#">NM_001353684</a> , <a href="#">NM_001353685</a> , <a href="#">NM_001353686</a> , <a href="#">NM_001353687</a> , <a href="#">NM_001353688</a> , <a href="#">NM_001353689</a> , <a href="#">NM_001353690</a> , <a href="#">NM_001353691</a> , <a href="#">NM_001353692</a> , <a href="#">NM_001353693</a> , <a href="#">NM_001353694</a> , <a href="#">NM_003253.1</a> , <a href="#">NM_003253.2</a> , <a href="#">BC117192</a> , <a href="#">BC117196</a> , <a href="#">BC143970</a> , <a href="#">BC143980</a> , <a href="#">NM_003253.3</a>
UniProt ID:	<a href="#">Q13009</a>



[View online »](#)

<b>Summary:</b>	<p>This gene encodes a RAC1-specific guanine nucleotide exchange factor (GEF). GEFs mediate the exchange of guanosine diphosphate (GDP) for guanosine triphosphate (GTP). The binding of GTP induces a conformational change in RAC1 that allows downstream effectors to bind and transduce a signal. This gene thus regulates RAC1 signaling pathways that affect cell shape, migration, adhesion, growth, survival, and polarity, as well as influencing actin cytoskeletal formation, endocytosis, and membrane trafficking. This gene thus plays an important role in cell invasion, metastasis, and carcinogenesis. In addition to RAC1, the encoded protein activates additional Rho-like GTPases such as CDC42, RAC2, RAC3 and RHOA. This gene encodes multiple protein isoforms that experience a diverse array of intramolecular, protein-protein, and phosphorylation interactions as well as phosphoinositide binding. Both the longer and shorter isoforms have C-terminal Dbl homology (DH) and pleckstrin homology (PH) domains while only the longer isoforms of this gene have the N-terminal myristoylation site and the downstream N-terminal PH domain, ras-binding domain (RBD), and PSD-95/DlgA/ZO-1 (PDZ) domain. [provided by RefSeq, Jul 2017]</p>
<b>shRNA Design:</b>	<p>These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a>.</p>
<b>Performance Guaranteed:</b>	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>