

Product datasheet for **TG311240**

NCoR (NCOR1) Human shRNA Plasmid Kit (Locus ID 9611)

Product data:

Product Type:	shRNA Plasmids
Product Name:	NCoR (NCOR1) Human shRNA Plasmid Kit (Locus ID 9611)
Locus ID:	9611
Synonyms:	hN-CoR; N-CoR; N-CoR1; PPP1R109; TRAC1
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	NCOR1 - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 9611). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	NM_001190438 , NM_001190440 , NM_006311 , NM_006311.1 , NM_006311.2 , NM_006311.3 , NM_001190440.1 , NM_001190438.1 , BC167431 , BC026028 , BC035748 , BC050594 , BC056862 , BC058511 , BC068996 , BC107772 , BC142648 , BC172437 , NM_006311.4
UniProt ID:	O75376
Summary:	This gene encodes a protein that mediates ligand-independent transcription repression of thyroid-hormone and retinoic-acid receptors by promoting chromatin condensation and preventing access of the transcription machinery. It is part of a complex which also includes histone deacetylases and transcriptional regulators similar to the yeast protein Sin3p. This gene is located between the Charcot-Marie-Tooth and Smith-Magenis syndrome critical regions on chromosome 17. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 17 and 20.[provided by RefSeq, Jun 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).