

## Product datasheet for **TG311699**

### LOXL2 Human shRNA Plasmid Kit (Locus ID 4017)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | shRNA Plasmids   |
| Product Name:             | LOXL2 Human shRNA Plasmid Kit (Locus ID 4017)  |
| Locus ID:                 | 4017   |
| Synonyms:                 | LOR; LOR2; WS9-14  |
| Vector:                   | pGFP-V-RS (TR30007)  |
| E. coli Selection:        | Kanamycin  |
| Mammalian Cell Selection: | Puromycin  |
| Format:                   | Retroviral plasmids  |
| Components:               | LOXL2 - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 4017).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.   |
| RefSeq:                   | <a href="#">NM_002318</a> , <a href="#">NM_002318.1</a> , <a href="#">NM_002318.2</a> , <a href="#">BC000594</a> , <a href="#">BC000594.2</a> , <a href="#">BC020481</a> , <a href="#">NM_002318.3</a>   |
| UniProt ID:               | <a href="#">Q9Y4K0</a>   |
| Summary:                  | This gene encodes a member of the lysyl oxidase gene family. The prototypic member of the family is essential to the biogenesis of connective tissue, encoding an extracellular copper-dependent amine oxidase that catalyses the first step in the formation of crosslinks in collagens and elastin. A highly conserved amino acid sequence at the C-terminus end appears to be sufficient for amine oxidase activity, suggesting that each family member may retain this function. The N-terminus is poorly conserved and may impart additional roles in developmental regulation, senescence, tumor suppression, cell growth control, and chemotaxis to each member of the family. [provided by RefSeq, Jul 2008] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .   |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).