

## Product datasheet for **TG312095**

### Integrin alpha 3 (ITGA3) Human shRNA Plasmid Kit (Locus ID 3675)

#### Product data:

|                           |   |
|---------------------------|---|
| Product Type:             | shRNA Plasmids  |
| Product Name:             | Integrin alpha 3 (ITGA3) Human shRNA Plasmid Kit (Locus ID 3675)  |
| Locus ID:                 | 3675  |
| Synonyms:                 | CD49C; FRP-2; GAP-B3; GAPB3; ILNEB; MSK18; VCA-2; VL3A; VLA3a   |
| Vector:                   | pGFP-V-RS (TR30007)   |
| E. coli Selection:        | Kanamycin   |
| Mammalian Cell Selection: | Puromycin   |
| Format:                   | Retroviral plasmids   |
| Components:               | ITGA3 - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 3675).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.  |
| RefSeq:                   | <a href="#">NM_002204</a> , <a href="#">NM_005501</a> , <a href="#">NM_002204.1</a> , <a href="#">NM_002204.2</a> , <a href="#">NM_002204.3</a> , <a href="#">NM_005501.1</a> ,<br><a href="#">NM_005501.2</a> , <a href="#">BC015344</a> , <a href="#">BC043917</a> , <a href="#">BC136636</a> , <a href="#">BC144328</a> , <a href="#">BC150190</a> , <a href="#">BM978115</a> ,<br><a href="#">NM_002204.4</a>   |
| UniProt ID:               | <a href="#">P26006</a>  |
| Summary:                  | The gene encodes a member of the integrin alpha chain family of proteins. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain that function as cell surface adhesion molecules. The encoded preproprotein is proteolytically processed to generate light and heavy chains that comprise the alpha 3 subunit. This subunit joins with a beta 1 subunit to form an integrin that interacts with extracellular matrix proteins including members of the laminin family. Expression of this gene may be correlated with breast cancer metastasis. [provided by RefSeq, Oct 2015] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .  |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).