

## Product datasheet for **TG313603**

### CYP24A1 Human shRNA Plasmid Kit (Locus ID 1591)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CYP24A1 Human shRNA Plasmid Kit (Locus ID 1591)
Locus ID:	1591
Synonyms:	CP24; CYP24; HCAI; HCINF1; P450-CC24
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CYP24A1 - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 1591). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	<a href="#">NM_000782</a> , <a href="#">NM_001128915</a> , <a href="#">NM_000782.1</a> , <a href="#">NM_000782.2</a> , <a href="#">NM_000782.3</a> , <a href="#">NM_000782.4</a> , <a href="#">NM_001128915.1</a> , <a href="#">BC109083</a> , <a href="#">BC109084</a> , <a href="#">BM928702</a> , <a href="#">NM_001128915.2</a> , <a href="#">NM_000782.5</a>
UniProt ID:	<a href="#">Q07973</a>
Summary:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This mitochondrial protein initiates the degradation of 1,25-dihydroxyvitamin D3, the physiologically active form of vitamin D3, by hydroxylation of the side chain. In regulating the level of vitamin D3, this enzyme plays a role in calcium homeostasis and the vitamin D endocrine system. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).