

## Product datasheet for **TG313961**

### Cofilin 1 (CFL1) Human shRNA Plasmid Kit (Locus ID 1072)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Cofilin 1 (CFL1) Human shRNA Plasmid Kit (Locus ID 1072)
Locus ID:	1072
Synonyms:	CFL; cofilin; HEL-S-15
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CFL1 - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 1072). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	<a href="#">NM_005507</a> , <a href="#">NM_005507.1</a> , <a href="#">NM_005507.2</a> , <a href="#">BC012318</a> , <a href="#">BC012318.1</a> , <a href="#">BC011005</a> , <a href="#">BC012265</a> , <a href="#">BC018256</a> , <a href="#">BC024975</a> , <a href="#">NM_005507.3</a>
UniProt ID:	<a href="#">P23528</a>
Summary:	The protein encoded by this gene can polymerize and depolymerize F-actin and G-actin in a pH-dependent manner. Increased phosphorylation of this protein by LIM kinase aids in Rho-induced reorganization of the actin cytoskeleton. Cofilin is a widely distributed intracellular actin-modulating protein that binds and depolymerizes filamentous F-actin and inhibits the polymerization of monomeric G-actin in a pH-dependent manner. It is involved in the translocation of actin-cofilin complex from cytoplasm to nucleus.[supplied by OMIM, Apr 2004]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).