

## Product datasheet for **TG314881**

### AHNAK Human shRNA Plasmid Kit (Locus ID 79026)

#### Product data:

|                           |   |
|---------------------------|---|
| Product Type:             | shRNA Plasmids  |
| Product Name:             | AHNAK Human shRNA Plasmid Kit (Locus ID 79026)  |
| Locus ID:                 | 79026   |
| Synonyms:                 | AHNAKRS; PM227  |
| Vector:                   | pGFP-V-RS (TR30007)   |
| E. coli Selection:        | Kanamycin   |
| Mammalian Cell Selection: | Puromycin   |
| Format:                   | Retroviral plasmids   |
| Components:               | AHNAK - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 79026).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.   |
| RefSeq:                   | <a href="#">NM_001620</a> , <a href="#">NM_024060</a> , <a href="#">NM_001346445</a> , <a href="#">NM_001346446</a> , <a href="#">NM_024060.1</a> , <a href="#">NM_024060.2</a> , <a href="#">NM_024060.3</a> , <a href="#">NM_001620.1</a> , <a href="#">BC000926</a> , <a href="#">BC012477</a> , <a href="#">BC017483</a> , <a href="#">BC128460</a> , <a href="#">NM_024060.4</a>   |
| UniProt ID:               | <a href="#">Q09666</a>  |
| Summary:                  | The protein encoded by this gene is a large (700 kDa) structural scaffold protein consisting of a central domain with 128 aa repeats. The encoded protein may play a role in such diverse processes as blood-brain barrier formation, cell structure and migration, cardiac calcium channel regulation, and tumor metastasis. A much shorter variant encoding a 17 kDa isoform exists for this gene, and the shorter isoform initiates a feedback loop that regulates alternative splicing of this gene. [provided by RefSeq, Oct 2016] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .  |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).