

## Product datasheet for **TG317170**

### ING4 Human shRNA Plasmid Kit (Locus ID 51147)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ING4 Human shRNA Plasmid Kit (Locus ID 51147)
Locus ID:	51147
Synonyms:	my036; p29ING4
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ING4 - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 51147). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	<a href="#">NM_001127582</a> , <a href="#">NM_001127583</a> , <a href="#">NM_001127584</a> , <a href="#">NM_001127585</a> , <a href="#">NM_001127586</a> , <a href="#">NM_016162</a> , <a href="#">NM_198287</a> , <a href="#">NM_016162.1</a> , <a href="#">NM_016162.2</a> , <a href="#">NM_016162.3</a> , <a href="#">NM_001127586.1</a> , <a href="#">NM_001127585.1</a> , <a href="#">NM_001127584.1</a> , <a href="#">NM_001127583.1</a> , <a href="#">NM_001127582.1</a> , <a href="#">BC007781</a> , <a href="#">BC007781.2</a> , <a href="#">BC013038</a> , <a href="#">BC095434</a> , <a href="#">BM727982</a> , <a href="#">BM922497</a> , <a href="#">NM_001127586.2</a> , <a href="#">NM_001127585.2</a> , <a href="#">NM_001127582.2</a> , <a href="#">NM_001127583.2</a> , <a href="#">NM_001127584.2</a>
UniProt ID:	<a href="#">Q9UNL4</a>
Summary:	This gene encodes a tumor suppressor protein that contains a PHD-finger, which is a common motif in proteins involved in chromatin remodeling. This protein can bind TP53 and EP300/p300, a component of the histone acetyl transferase complex, suggesting its involvement in the TP53-dependent regulatory pathway. Multiple alternatively spliced transcript variants have been observed that encode distinct proteins. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).