

## Product datasheet for **TG321436**

### CD24 Human shRNA Plasmid Kit (Locus ID 100133941)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CD24 Human shRNA Plasmid Kit (Locus ID 100133941)
Locus ID:	100133941
Synonyms:	CD24A
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CD24 - Human, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 100133941). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	<a href="#">BC064619</a> , <a href="#">NM_001291737</a> , <a href="#">NM_001291738</a> , <a href="#">NM_001291739</a> , <a href="#">NM_013230</a> , <a href="#">NR_117089</a> , <a href="#">NR_117090</a> , <a href="#">NM_001359084</a> , <a href="#">NM_013230.2</a> , <a href="#">NM_013230.3</a> , <a href="#">NM_001291737.1</a> , <a href="#">NM_001291738.1</a> , <a href="#">NM_001291739.1</a> , <a href="#">BC064619.1</a>
UniProt ID:	<a href="#">P25063</a>
Summary:	This gene encodes a sialoglycoprotein that is expressed on mature granulocytes and B cells and modulates growth and differentiation signals to these cells. The precursor protein is cleaved to a short 32 amino acid mature peptide which is anchored via a glycosyl phosphatidylinositol (GPI) link to the cell surface. This gene was missing from previous genome assemblies, but is properly located on chromosome 6. Non-transcribed pseudogenes have been designated on chromosomes 1, 15, 20, and Y. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).