

## Product datasheet for **TG500064**

### Ahr Mouse shRNA Plasmid (Locus ID 11622)

#### Product data:

|                           |   |
|---------------------------|---|
| Product Type:             | shRNA Plasmids  |
| Product Name:             | Ahr Mouse shRNA Plasmid (Locus ID 11622)  |
| Locus ID:                 | 11622   |
| Synonyms:                 | A; Ah; Ahh; Ahre; bHLHe7; bHLHe76; In   |
| Vector:                   | pGFP-V-RS (TR30007)   |
| E. coli Selection:        | Kanamycin   |
| Mammalian Cell Selection: | Puromycin   |
| Format:                   | Retroviral plasmids   |
| Components:               | Ahr - Mouse, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 11622). 5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.  |
| RefSeq:                   | <a href="#">NM_013464</a> , <a href="#">NM_013464.1</a> , <a href="#">NM_013464.2</a> , <a href="#">NM_013464.3</a> , <a href="#">NM_013464.4</a> , <a href="#">BC040248</a> , <a href="#">BC152973</a> , <a href="#">BC156629</a>  |
| UniProt ID:               | <a href="#">P30561</a>  |
| Summary:                  | The protein encoded by this gene is a ligand-activated helix-loop-helix transcription factor involved in the regulation of biological responses to planar aromatic hydrocarbons. This receptor has been shown to regulate xenobiotic-metabolizing enzymes such as cytochrome P450. Before ligand binding, the encoded protein is sequestered in the cytoplasm; upon ligand binding, this protein moves to the nucleus and stimulates transcription of target genes. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2015] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .  |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).