

## Product datasheet for **TG500130**

### Arf6 Mouse shRNA Plasmid (Locus ID 11845)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Arf6 Mouse shRNA Plasmid (Locus ID 11845)
Locus ID:	11845
Synonyms:	AI788669; AW496366
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Arf6 - Mouse, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 11845). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	<a href="#">BC003478</a> , <a href="#">BC083112</a> , <a href="#">NM_007481</a> , <a href="#">NM_007481.1</a> , <a href="#">NM_007481.2</a> , <a href="#">NM_007481.3</a>
UniProt ID:	<a href="#">P62331</a>
Summary:	GTP-binding protein involved in protein trafficking that regulates endocytic recycling and cytoskeleton remodeling (PubMed:11950392). Required for normal completion of mitotic cytokinesis. Involved in the regulation of dendritic spine development, contributing to the regulation of dendritic branching and filopodia extension. Plays an important role in membrane trafficking, during junctional remodeling and epithelial polarization. Regulates surface levels of adherens junction proteins such as CDH1 (PubMed:29420262, PubMed:20080746). Required for NTRK1 sorting to the recycling pathway from early endosomes (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).