

Product datasheet for **TG501475**

Nin Mouse shRNA Plasmid (Locus ID 18080)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Nin Mouse shRNA Plasmid (Locus ID 18080)
Locus ID:	18080
Synonyms:	3110068G20Rik; AI385615; AU024711; mKIAA1565
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Nin - Mouse, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 18080). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	NM_001081453 , NM_001286079 , NM_001286080 , NM_008697 , NR_104397 , NR_104398 , NM_001081453.1 , NM_008697.1 , NM_008697.2 , NM_008697.3 , NM_008697.4 , NM_001286080.1 , NM_001286080.2 , NM_001286079.1 , NM_001286079.2 , BC144913 , BC137789 , BC144914
UniProt ID:	Q61043
Summary:	Centrosomal protein required for the positioning and anchorage of the microtubule minus-end in epithelial cells (PubMed:15784680, PubMed:10934040). May also act as a centrosome maturation factor (By similarity). May play a role in microtubule nucleation, by recruiting the gamma-tubulin ring complex to the centrosome (PubMed:15784680). Overexpression does not perturb nucleation or elongation of microtubules but suppresses release of microtubules (By similarity). Required for centriole organization and microtubule anchoring at the mother centriole (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).