

Product datasheet for **TG509852**

Hnrnpk Mouse shRNA Plasmid (Locus ID 15387)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Hnrnpk Mouse shRNA Plasmid (Locus ID 15387)
Locus ID:	15387
Synonyms:	hn; Hnrpk; KBB; KBBP; NOVA
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Hnrnpk - Mouse, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 15387). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	BC006694 , BC089328 , NM_001301340 , NM_001301341 , NM_001301343 , NM_001301344 , NM_001301345 , NM_025279 , NR_153501 , NR_153502 , NR_153503 , NR_153504 , NM_025279.1 , NM_025279.2 , NM_025279.3 , NM_001301345.1 , NM_001301344.1 , NM_001301340.1 , NM_001301341.1 , NM_001301343.1
UniProt ID:	P61979
Summary:	The protein encoded by this gene is a component of the heterogeneous nuclear ribonucleoprotein (hnRNP) complex, and is a poly-cytosine binding protein (PCBP). It is found in multiple subcellular compartments including the nucleus, cytoplasm and mitochondria. This gene product is thought to interact with RNA, DNA, and protein, and is involved in multiple cellular processes, including transcription, chromatin remodeling, DNA damage response, signal transduction, mRNA splicing, export, and translation. Multiple transcript variants and protein isoforms exist, with some isoforms containing a unique C-terminus. There are four pseudogenes of this gene, found on chromosomes 2, 3, 7 and 13. [provided by RefSeq, Aug 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).