

## Product datasheet for **TG513177**

### Zeb1 Mouse shRNA Plasmid (Locus ID 21417)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | shRNA Plasmids   |
| Product Name:             | Zeb1 Mouse shRNA Plasmid (Locus ID 21417)  |
| Locus ID:                 | 21417  |
| Synonyms:                 | 3110032K11Rik; AREB6; BZP; MEB1; Nil2; TCF-8; Tcf8; Tcf18; Tw; ZEB; Zfh1a; Zfx1a; Zfx1ha   |
| Vector:                   | pGFP-V-RS (TR30007)  |
| E. coli Selection:        | Kanamycin  |
| Mammalian Cell Selection: | Puromycin  |
| Format:                   | Retroviral plasmids  |
| Components:               | Zeb1 - Mouse, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 21417).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.   |
| RefSeq:                   | <a href="#">BC139768</a> , <a href="#">BC139769</a> , <a href="#">NM_011546</a> , <a href="#">NM_001360981</a> , <a href="#">NM_001360982</a> , <a href="#">NM_011546.1</a> , <a href="#">NM_011546.2</a> , <a href="#">NM_011546.3</a> , <a href="#">BC090987</a>   |
| UniProt ID:               | <a href="#">Q64318</a>   |
| Summary:                  | Acts as a transcriptional repressor. Binds to E-box sequences in the immunoglobulin heavy chain enhancer as well as in the regulatory regions of many other tissue-specific genes. Represses E-cadherin promoter and induces an epithelial-mesenchymal transition (EMT) by recruiting SMARCA4/BRG1. Represses BCL6 transcription in the presence of the corepressor CTBP1 (By similarity). Positively regulates neuronal differentiation. Represses RCOR1 transcription activation during neurogenesis. Represses transcription by binding to the E box (5'-CANNTG-3'). Promotes tumorigenicity by repressing stemness-inhibiting microRNAs. [UniProtKB/Swiss-Prot Function] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .   |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).