

Product datasheet for **TG513741**

Twist1 Mouse shRNA Plasmid (Locus ID 22160)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Twist1 Mouse shRNA Plasmid (Locus ID 22160)
Locus ID:	22160
Synonyms:	bHLHa; bHLHa38; M-Twi; M-Twist; pd; Pde; pdt; Pluri; Ska; Ska10; Ską Twist
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Twist1 - Mouse, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 22160). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	BC033434 , BC083139 , NM_011658 , NM_011658.1 , NM_011658.2
UniProt ID:	P26687
Summary:	Basic helix-loop-helix (bHLH) transcription factors have been implicated in cell lineage determination and differentiation. This gene encodes a bHLH transcription factor that is evolutionarily conserved from invertebrates to humans, and was originally identified in <i>Drosophila</i> as an essential gene involved in early mesoderm development and dorsal-ventral patterning in the embryo. This protein plays a role in cancer by regulating the epithelial-mesenchymal transition (EMT), a process that is critical for metastasis initiation, and promoting tumor progression. Mutations in the human gene are associated with Saethre-Chotzen syndrome (SCS). Mice with heterozygous mutations in this gene exhibit cranofacial and structural defects similar to those seen in human SCS patients. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).