

Product datasheet for **TG514066**

Tcf7 Mouse shRNA Plasmid (Locus ID 21414)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Tcf7 Mouse shRNA Plasmid (Locus ID 21414)
Locus ID:	21414
Synonyms:	AI465550; TCF-1; Tcf1
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Tcf7 - Mouse, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 21414). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	NM_009331 , NM_009331.1 , NM_009331.2 , NM_009331.3 , NM_009331.4 , BC138596 , BC056630 , BC145343
UniProt ID:	Q00417
Summary:	This gene encodes a transcription factor which is a member of the T-cell specific transcription factor family. The encoded protein is distinct from the hepatic transcription factor, transcription factor 1, which is also referred to by the symbol Tcf1. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).