

## Product datasheet for **TG516934**

### Nmnat1 Mouse shRNA Plasmid (Locus ID 66454)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Nmnat1 Mouse shRNA Plasmid (Locus ID 66454)
Locus ID:	66454
Synonyms:	2610529L11Rik; 5730441G13Rik; D4Cole1e; nmnat
Vector:	pGFP-V-RS (TR30007)
E. coli Selection:	Kanamycin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Nmnat1 - Mouse, 4 unique 29mer shRNA constructs in retroviral GFP vector(Gene ID = 66454). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-V-RS Vector, TR30013, included for free.
RefSeq:	<a href="#">NM_133435</a> , <a href="#">NM_001356357</a> , <a href="#">NM_133435.1</a> , <a href="#">BC138500</a> , <a href="#">BC138502</a> , <a href="#">NM_133435.2</a>
UniProt ID:	<a href="#">Q9EPA7</a>
Summary:	Catalyzes the formation of NAD(+) from nicotinamide mononucleotide (NMN) and ATP (PubMed:15381699). Can also use the deamidated form; nicotinic acid mononucleotide (NaMN) as substrate with the same efficiency (By similarity). Can use triazofurin monophosphate (TrMP) as substrate (By similarity). Also catalyzes the reverse reaction, i.e. the pyrophosphorolytic cleavage of NAD(+) (By similarity). For the pyrophosphorolytic activity, prefers NAD(+) and NaAD as substrates and degrades NADH, nicotinic acid adenine dinucleotide phosphate (NAD) and nicotinamide guanine dinucleotide (NGD) less effectively (By similarity). Involved in the synthesis of ATP in the nucleus, together with PARP1, PARG and NUDT5 (By similarity). Nuclear ATP generation is required for extensive chromatin remodeling events that are energy-consuming (By similarity). Fails to cleave phosphorylated dinucleotides NADP(+), NADPH and NaADP(+) (By similarity). Protects against axonal degeneration following mechanical or toxic insults (PubMed:15310905, PubMed:16914673). Delays axonal degeneration after axotomy. Results in a >10-fold increase in intact neurites 72 hours after injury (PubMed:16914673).[UniProtKB/Swiss-Prot Function]



[View online »](#)

<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .
<b>Performance Guaranteed:</b>	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at <a href="mailto:techsupport@origene.com">techsupport@origene.com</a>. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>