

Product datasheet for **TL300058**

ZNF683 Human shRNA Plasmid Kit (Locus ID 257101)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ZNF683 Human shRNA Plasmid Kit (Locus ID 257101)
Locus ID:	257101
Synonyms:	Hobit
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ZNF683 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 257101). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001114759 , NM_001307925 , NM_173574 , NM_173574.1 , NM_173574.3 , NM_001114759.1 , NM_001114759.2 , BC029505 , BC029505.1 , BC028731 , NM_001114759.3
UniProt ID:	Q8IZ20
Summary:	Transcription factor that mediates a transcriptional program in various innate and adaptive immune tissue-resident lymphocyte T-cell types such as tissue-resident memory T (Trm), natural killer (trNK) and natural killer T (NKT) cells and negatively regulates gene expression of proteins that promote the egress of tissue-resident T-cell populations from non-lymphoid organs. Plays a role in the development, retention and long-term establishment of adaptive and innate tissue-resident lymphocyte T cell types in non-lymphoid organs, such as the skin and gut, but also in other nonbarrier tissues like liver and kidney, and therefore may provide immediate immunological protection against reactivating infections or viral reinfection. Plays also a role in the differentiation of both thymic and peripheral NKT cells. Negatively regulates the accumulation of interferon-gamma (IFN-gamma) in NKT cells at steady state or after antigenic stimulation. Positively regulates granzyme B production in NKT cells after innate stimulation. Associates with the transcriptional repressor PRDM1/BLIMP1 to chromatin at gene promoter regions.[UniProtKB/Swiss-Prot Function]


[View online »](#)

shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>