

## Product datasheet for **TL300438**

### TAZ (WWTR1) Human shRNA Plasmid Kit (Locus ID 25937)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	TAZ (WWTR1) Human shRNA Plasmid Kit (Locus ID 25937)
Locus ID:	25937
Synonyms:	TAZ
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	WWTR1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 25937). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001168278</a> , <a href="#">NM_001168280</a> , <a href="#">NM_015472</a> , <a href="#">NM_001348362</a> , <a href="#">NM_015472.1</a> , <a href="#">NM_015472.2</a> , <a href="#">NM_015472.3</a> , <a href="#">NM_015472.4</a> , <a href="#">NM_015472.5</a> , <a href="#">NM_001168278.1</a> , <a href="#">NM_001168278.2</a> , <a href="#">NM_001168280.1</a> , <a href="#">NM_001168280.2</a> , <a href="#">BC014052</a> , <a href="#">BC014052.2</a> , <a href="#">NM_001168280.3</a> , <a href="#">NM_015472.6</a>
UniProt ID:	<a href="#">Q9GZV5</a>
Summary:	Transcriptional coactivator which acts as a downstream regulatory target in the Hippo signaling pathway that plays a pivotal role in organ size control and tumor suppression by restricting proliferation and promoting apoptosis. The core of this pathway is composed of a kinase cascade wherein STK3/MST2 and STK4/MST1, in complex with its regulatory protein SAV1, phosphorylates and activates LATS1/2 in complex with its regulatory protein MOB1, which in turn phosphorylates and inactivates YAP1 oncoprotein and WWTR1/TAZ. WWTR1 enhances PAX8 and NKX2-1/TTF1-dependent gene activation. Regulates the nuclear accumulation of SMADS and has a key role in coupling them to the transcriptional machinery such as the mediator complex. Regulates embryonic stem-cell self-renewal, promotes cell proliferation and epithelial-mesenchymal transition.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).